

Distillation & Walk-in Fume Cupboards

Expelling all
doubts



Introduction

The Clean Air Distillation and Walk-in fume cupboards are ideal for use with distillation rigs, tall equipment or where full access is required.

Clean Air Distillation & Walk-in range

Available in the standard widths of 1200mm, 1500mm, 1800mm and 2000mm wide, with larger or smaller units available on request.

Aerodynamically designed.

Choice of chamber construction.

Choice of low-level and intermediate worktops, which are fully retractable.

Independent twin-sash arrangement.

Distillation standard chamber height of 1555mm (BTH @ 400mm.)

Walk-in chamber height of 1955mm.

High, low and intermediate level baffles to purge both lighter and heavier than air fumes.

Performance and safety

The Clean Air Distillation and Walk-in fume cupboards are certified to BS7258.



Technical & Design

Twin-skin shell in distinctive styling and a choice of colours and finishes.

Sprung PVC guide channels for directional stability, with smooth rise and fall sash action.

Twin sashes counterbalanced using double stainless steel cables and ball race pulleys of the fail-safe design.

Adaptable for use with energy saving Variable Air Volume (VAV) systems.

Matching cosmetic infill panels from top of unit to ceiling to shroud ducting.

Contact

Clean Air Limited
Dunscar Industrial Estate
Blackburn Road
Bolton
Lancashire
BL7 9PQ
United Kingdom
Tel +44 (0) 1204 591115
Fax +44 (0) 1204 591116
Website: cleanairltd.co.uk
E-Mail: sales@cleanairltd.co.uk

2000mm wide Walk-in unit shown with retractable worktop.

Features

- Choice of retractable worktop heights from 200mm to 400mm.



- Drip cups & services available.



Annual Maintenance

Clean Air offer annual maintenance and inspection of both new and existing fume cupboards and extract systems.

