

2E. Education Fixed Fume cupboards



Standard 1000mm wide unit shown with infill panels and steel cupboards.

Features

- Clean and attractive design.
- 16swg Zintec epoxy powder coated construction White Ral 9010 textured finish.
- High visibility and ideal for demonstrations.
- Available in Fixed or Semi-mobile design.
- Fan switch with indicator.
- Airflow control panel.
- Unique sash lock at 400mm.
- Light switch with indicator.
- Double switched socket outlet.
- Water outlet – remote operated.
- Gas outlet – remote operated.
- Vulcathene circular drip cup.

Introduction

The Clean Air Educational range is a high quality, robust design of fume cupboard at a competitive price.

Clean Air Education range

Available in the standard widths of 1000mm, 1200mm & 1500mm.

Type tested to BS7258 with unique ultimate levels of containment.

Aerodynamically designed.

Mechanical water and gas services are controlled with remote valves and outlets are mounted off the worktop to prevent corrosion from spillages.

Unique sash locking device at 400mm for operator protection.

Air baffles to ensure optimum containment levels on site.

Sash handle profile is designed to improve containment within the work chamber and eliminate eddies.

Dished seamless high quality GRP worktop.

Steel lockable storage unit included.

ducted,
filtration
and
contamination
control
equipment

Performance and safety

The Clean Air **Education** range of Fume Cupboards have been independently tested to: -

CLEAPSS

BS7258

Building Bulletin 88 (BB88)

Design Note 29 (DN29)

SSERC

Technical & Design

Full height toughened glass sides and sash for demonstrations.

All the safety features and aerodynamics of an Industrial design of fume cupboard.

Matching cosmetic infill panels from top of unit to ceiling to shroud ducting.

Warranty & after sales

Clean Air offer a full one-year UNCONDITIONAL warranty on all parts and labour.

Clean Air Ltd. Unit 16 Dunsar
Industrial Estate, Blackburn
Road, Bolton. BL7 9PQ.
Tel: +44 (0)1204 591115.
Fax: +44 (0)1204 591116.

