

# 1E. Education SSR2000 Mobile Filtration Fume cupboard



Standard 1000mm wide unit shown with water and gas services.

## Features

- Clean and attractive design.
- Super rigid, lightweight aluminium construction.
- 360 degree visibility.
- Operator protection by fail safe toughened glass sliding sash.
- No harmful fumes to atmosphere – integral 18kg carbon filter.
- Airflow control panel membrane.
- Unique sash lock at 400mm.
- Light switch with indicator.
- RCD Protected double socket.
- Water outlet – remote operated.
- Gas outlet – remote operated.
- Drip cup, bottle trap & 2m hose.
- Protective pre-filter.

## Introduction

The **SSR2000** is a fully mobile Filtration Fume cupboard, ideal for Education & demonstrations.

## Sliding Sash Recirculating Fume Cupboard

Quick release connections for instant installation and easy relocation.

Twin lights give 1000 Lux, with 360-degree visibility for demonstrations

Minimum airflow of 0.4m/s at a sash height of 400mm.

Designed to pass through a standard door opening.

Unique sash locking device at 400mm for operator protection.



Recessed controls.

Very low noise level – only 50dBA.

Dished seamless high quality GRP worktop.

Ideal for listed buildings or where ductwork is difficult to install.

Full height toughened glass sides and sash for demonstrations.

ducted,  
filtration  
and  
contamination  
control  
equipment

## Performance and safety

The Clean Air **SSR2000** Fume Cupboard has been independently tested to : -

**CLEAPSS**

**BS7258**

**Building Bulletin 88 (BB88)**

**Design Note 29 (DN29)**

**SSERC**

## Technical & Design



Safe & easy filter change.

1000mm wide x 1950mm high x 700mm deep.

Energy conservation.

## Warranty & after sales

Clean Air offer a full one-year UNCONDITIONAL warranty on all parts and labour.

Clean Air Ltd. Unit 16 Dunscair Industrial Estate, Blackburn Road, Bolton. BL7 9PQ.  
Tel: +44 (0)1204 591115.  
Fax: +44 (0)1204 591116.

