

# 3E. Education Semi-mobile Fume cupboards



1000mm wide unit shown with steel cupboards, airflow control panel and flexible / rigid ductwork connection.

## Features

- Clean and attractive design.
- 16swg Zintec epoxy powder coated construction White Ral 9010 textured finish.
- Glazed to 360-degrees, ideal for demonstrations.
- Airflow control panel.
- Unique sash lock at 400mm.
- Light switch with indicator.
- Swivel castors with the front two being lockable.
- Double switched socket outlet.
- Water outlet – remote operated.
- Gas outlet – remote operated.
- Vulcathene circular drip cup.

## Introduction

The Education range is a high quality, robust design of fume cupboard at a competitive price.

## Education range

Available in the standard widths of 1000mm, 1200mm & 1500mm.

Type tested to BS7258 with unique ultimate levels of containment.

Aerodynamically designed.

Mechanical water and gas services are controlled with remote valves and outlets are mounted off the worktop to prevent corrosion from spillages.

Unique sash locking device at 400mm for operator protection.

Air baffles to ensure optimum containment levels on site.

Sash handle profile is designed to improve containment within the work chamber and eliminate eddies.

Dished seamless high quality GRP worktop.

Steel lockable storage unit included.

Complete with 2m long water, gas and waste hoses.

ducted,  
filtration  
and  
contamination  
control  
equipment

## Performance and safety

The Education range of Fume Cupboards have been independently tested to: -

**CLEAPSS**

**BS7258**

**Building Bulletin 88 (BB88)**

**Design Note 29 (DN29)**

**SSERC**

## Technical & Design

Full height toughened glass sides, back, rear baffles and sash for demonstrations.

All the safety features and aerodynamics of an Industrial design of fume cupboard.

Flexible ductwork to allow the unit to be pulled out 1.5m.

## Warranty & after sales

All fume cupboards include a full one-year UNCONDITIONAL warranty on all parts and labour.

Clean Air Ltd. Unit 16 Dunscair  
Industrial Estate, Blackburn  
Road, Bolton. BL7 9PQ.  
Tel: +44 (0)1204 591115.  
Fax: +44 (0)1204 591116.

