

High Sash Fume Cupboard

Expelling all
doubts



2000mm wide unit shown with trough services and vented steel cupboards.

Features

- Larger front viewing height and sash opening



- Increased internal chamber heights for taller equipment.
- Increased height for maximum quantities of services



- Scaffold grids, scaffold bosses and tissue guards available.

Introduction

The Clean Air High Sash fume cupboard has been designed to accommodate maximum quantities of services, together with a taller internal height and increased front viewing height.

Clean Air High Sash range

Available in the standard widths of 1200mm, 1500mm, 1800mm and 2000mm wide, with larger or smaller units available on request.

Aerodynamically designed.

Interior chamber smooth and easy to clean and maintain.

Unique sash locking device for operator protection.

Adjustable baffles to ensure optimum containment levels on site.

Choice of lining materials.

Choice of worktop materials.

6mm thick toughened glass sash.

Factory pre-plumbed and pre-wired for simple site connection.

Available in either bench mounted, Distillation or Walk-in version.

Available in either ducted or filtration format.

Annual Maintenance

Clean Air offer annual maintenance and inspection of both new and existing fume cupboards and extract systems.

Performance and safety

The Clean Air High Sash range of fume cupboards are certified to BS7258.

All Clean Air equipment conforms to the latest electrical safety requirements & EC directive for EMC.

Technical & Design

Twin-skin shell in distinctive styling and a choice of colours and finishes.

Sprung PVC guide channels for directional stability, with smooth rise and fall sash action.

Sash counterbalance using double stainless steel cables and ball race pulleys of the fail-safe design.

Adaptable for use with energy saving Variable Air Volume (VAV) systems.

Matching cosmetic infill panels from top of unit to ceiling to shroud ducting.

Contact

Clean Air Limited
Dunscar Industrial Estate
Blackburn Road
Bolton
Lancashire
BL7 9PQ
United Kingdom
Tel +44 (0) 1204 591115
Fax +44 (0) 1204 591116
Website: cleanairltd.co.uk
E-Mail: sales@cleanairltd.co.uk

