

# Standard Fume Cupboard



1500mm wide unit shown with infill panels and vented steel cupboards.

## Features

- Services are side wall mounted



- Choice of Mechanical & Electrical services



- Choice of chamber materials



- Sinks and drip cups available in a full range of sizes.

## Introduction

The Clean Air Standard fume cupboard is a high quality, robust design of fume cupboard at a competitive price, available in either ducted or filtration format.

## Clean Air Standard Range

Available in the standard widths of 1200mm, 1500mm, 1800mm and 2000mm wide, with larger or smaller units available on request.

Aerodynamically designed.

Interior chamber smooth and easy to clean and maintain.

Unique sash locking device for operator protection.

Adjustable baffles to ensure optimum containment levels on site.

Sash handle profile is designed to improve containment within the work chamber and eliminate eddies.

Factory pre-plumbed and pre-wired for rapid site connection.

Available in either bench mounted, Distillation or Walk-in version.

6mm thick toughened glass sash.

Choice of worktop materials.

Choice of lining materials.

Sash counterbalance using double stainless steel cables and ball race pulleys of the fail-safe design.

Matching cosmetic infill panels from top of unit to ceiling to shroud ducting.

## Annual Maintenance

Clean Air offer annual maintenance and inspection of both new and existing fume cupboards and extract systems.

Expelling all doubts



## Performance & Safety

The Clean Air range of Standard Fume Cupboards are certified to BS EN 14175

## Technical & Design

Twin-skin shell in distinctive styling and a choice of colours and finishes.

Sprung PVC guide channels for directional stability, with smooth rise and fall sash action.



## Contact

Clean Air Limited  
Dunscar Industrial Estate  
Blackburn Road  
Bolton  
Lancashire  
BL7 9PQ  
United Kingdom  
Tel +44 (0) 1204 591115  
Fax +44 (0) 1204 591116  
Website: [cleanairltd.co.uk](http://cleanairltd.co.uk)  
E-Mail: [sales@cleanairltd.co.uk](mailto:sales@cleanairltd.co.uk)

ducted,  
filtration  
and  
contamination  
control  
equipment



Certified to ISO 9001:2000  
Certificate Number GB006900